

### The Technical Information of Tablet-type Oxygen Indicators

Weight of tablet oxygen indicators (gram)	0.19±0.02	
Color-change Status	More than 0.5% oxygen concentration(aerobic conditions)	It shows blue
	Less than 0.1 % oxygen concentration ( Anaerobic conditions )	It shows red
	oxygen concentration: 0.5~0.1%	It shows Color transition
The time from blue to red	Less than 0.1 % oxygen concentration	2 ~ 3 hours
The time from red to blue	From Anaerobic conditions to aerobic conditions	3~15 minutes
The frequency of reversible color change	Not less than 4 times	
The temperature and humidity range for usage	2°C~40°C	Relative Humidity 30%~95%
Storage conditions	Low temperature, Anaerobic, far away from sunlight	
Safety	Nontoxic	